GERY

1A200-HD

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EN Instruction for use

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NO Brukerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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GERY - 1A200-HD

Instructions for use for Prosthetists/Orthotists
Please read before using or working on the device

1A20099
2021-07

Pass on § 3, 7, 8 and 9 of these instructions to the Patient

1. COMPONENTS INCLUDED

Designation	Ref.	Included/Sold separately
Foot shell	1A200-xx	Included
Malleolar clip	1D101	Included

2. DESCRIPTION, PROPERTIES AND MECHANISM OF ACTION

A. Description

A prosthetic foot designed to allow older adult patients to perform their everyday activities, consisting of:

- A core and foot shell which cannot be disassembled
- A malleolar clip

The hole in the lower face of the foot is used to identify the location of the gap between the big toe and the other toes, so that a separation can be made if required to allow sandal-type footwear to be worn.

To make the cut, use a band saw to saw through to the hole and open up the cut as required.

B. Properties

Weight (Size 25)	395 g	
Side	Right or Left	
Size	22 to 25	26 to 29
Malleolar clip	1D101-2225	1D101-2629
Maximum patient weight (Including carried load)	150 kg	

3. INTENDED USERS/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists/orthotists) who then instruct the patient in how to use it. The prescription is drawn up by a doctor. The doctor works with a prosthetist/orthotist to assess whether the patient is suited to using the device.

⚠ This device should be used on **ONE PATIENT ONLY**. It must not be reused on other patients.

This device is exclusively intended to be used as orthopaedic equipment for a below-knee or above-knee amputee whose activity level is low (1) to medium (2).

Maximum weight (including carried load): 150 kg
 Designed to accommodate a heel height of 0 to 20 mm.

⚠ **CAUTION: the GERY foot is not suitable for bilateral use by double amputees**

4. CLINICAL BENEFITS

The device provides optimal comfort for below-knee and above-knee amputees when they perform their everyday activities.


5. ACCESSORIES AND COMPATIBILITIES

The foot incorporates a male pyramid link that allows it to be used with female pyramid connectors (refer to our catalogue). The heel height can be adjusted from 0 to 20 mm.

6. ASSEMBLY AND FITTING TO THE PATIENT

A. Assembly and Aligning

Assemble in the normal way, checking that the load line passes 5 to 15 mm in front of the axis of the pyramid and provides a stable upright position for the patient.

The dynamic alignment is performed in accordance with good professional practices.

B. Finishing

Comply with traditional good professional practices in terms of fitting the malleolar clip to cosmetic elements.

7. DETECTING ANY MALFUNCTIONING

⚠ If you notice that the device is behaving abnormally or if its characteristics feel different, or if it has been badly knocked, consult your prosthetist/orthotist.

8. WARNINGS, CONTRAINDICATIONS, SIDE EFFECTS

A. Warnings

⚠ Using the foot in a way that goes against the advice given by your prosthetist/orthotist may damage or degrade the foot's components (e.g. carrying heavy loads, excessive stresses, use after its service life has expired, etc.)

⚠ **The device is only resistant to fresh water and to bad weather conditions, and requires compliance with the care instructions specified in §9.A.**

B. Contraindications

⚠ This foot is not intended to be used for activities associated with a risk of significant impact or excessive loading.

⚠ This foot is not resistant to salt water or chlorinated water.

C. Side effects

There are no side effects directly associated with the device.

Any serious incident that occurs which relates to the device must be reported to the manufacturer and to the competent authority of the member state in which the incident occurred.

9. MAINTENANCE, STORAGE, DISPOSAL AND SERVICE LIFE

A. Maintenance/cleaning

Clean the foot shell with a damp cloth or sponge.

⚠ All work on the foot should be performed by the prosthetist/orthotist; no other operations, such as lubrication, work on the fastenings or other work, need to be performed.

⚠ After being used in fresh water or in bad weather conditions

- Rinse the foot with clean water
- Dry thoroughly

The malleolar clip may need to be replaced (in the event of wear, loss, etc.)

B. Storage

Service and storage temperature: -20°C to +60°C

Relative humidity of the air: no restrictions

C. Disposal

The various components of the foot are considered as special waste: stainless steel, thermoplastic (polyamide, polyethylene), polyurethane. They must be treated in accordance with current, applicable legislation.

D. Service life

It is advisable to arrange for a prosthetist/orthotist to check the product once a year.

10. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of first declaration of conformity
--	--------------	--	-----------------	--	--

11. REGULATORY INFORMATION

This product is a CE-marked medical device that is certified as conforming with regulation (EU) 2017/745

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PROTEOR SAS