RUSH ROVER®

ROV

RUSH Kid

KID

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EN Instruction for use

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NO Bruerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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RUSH ROVER® , RUSH Kid
Instructions for use for prosthetists
Read before use

IFU-02-013
 Rev. A
 2023-08

Pass on § 10 to 17 of these instructions to the patient.

1. INCLUDED ITEMS

Part description	Part number	Included / Sold separately
RUSH ROVER	ROV-xx-x-RU*	Included
RUSH Kid	KID-xx-xP-RU*	Included
Spectra sock	S0-NPS-200xx-00*	Suitable sock included
Foot shell	FSK	Sold separately

* see catalog

2. DESCRIPTION

RUSH ROVER and Kid are energy-return prosthetic feet consisting of:

- A low-profile fiberglass composite sole plate
- A male pyramid connection
- A Spectra sock

3. PROPERTIES

Model	RUSH ROVER			RUSH Kid
Weight*	950 g / 2.09 lbs			505 g / 1.11 lbs
Size	22 to 24 cm	25 to 27 cm	28 to 30 cm	19 to 21 cm
Build height*	54 mm / 2.13"	60 mm / 2.38"	64 mm / 2.50"	51 mm / 2.0"
Heel height	10 mm / 3/8"			

*Weight based on size 26 (ROVER) and size 20 (Kid), cat. 4, with foot shell and Spectra sock

Build height based on sizes 23, 26 or 29 (ROVER) and size 20 (Kid), cat. 4, with foot shell, Spectra sock, and 10 mm heel height

RUSH ROVER has been tested according to ISO 10328 for a maximum patient weight up to 166 kg for 2 million cycles.

Selection of foot category based on patient's weight and impact level – RUSH ROVER										
Weight ^{*)}	lbs	0-105	106-140	141-175	176-210	211-245	246-280	281-315	316-350	351-365
	kg	0-48	49-64	65-79	80-95	96-111	112-127	128-143	144-159	160-166
Impact level	Low	1	2	3	4	5	6	7	8	9
	Moderate	1	2	3	4	5	6	7	8	9
	High	2	3	4	5	6	7	8	9	-

^{*)} Body mass limit not to be exceeded (ISO 10328)

RUSH Kid has been tested according to ISO 10328 for a maximum patient weight up to 61 kg for 2 million cycles.

Selection of foot category based on patient's weight and impact level – RUSH Kid										
Weight ^{*)}	lbs	0-30	31-45	46-60	61-75	76-90	91-105	106-120	121-135	
	kg	0-14	15-20	21-27	28-34	35-41	42-48	49-54	55-61	
Impact level	Low	1P	2P	3P	4P	5P	6P	7P	8P	
	Moderate	1P	2P	3P	4P	5P	6P	7P	8P	
	High	2P	3P	4P	5P	6P	7P	8P	9P	

^{*)} Body mass limit not to be exceeded (ISO 10328)

4. MECHANISM OF ACTION


The half rocker sole plate provides constant contact to the ground while eliminating dead spots, providing a smooth rollover, and natural energy return. The composite keel and sole plate store energy during initial contact and loading response which is then released during midstance. The keel and sole plate absorb energy during terminal stance and release the stored energy during pre-swing.

5. INTENDED USE/INDICATIONS

These medical devices are supplied to healthcare professionals (prosthetists), who will train the patient in their use. The prescription is drawn up by a doctor who assesses the patient's ability to use the device.

⚠ These devices are for multiple use on a **SINGLE PATIENT**. They must not be used on another patient.

These devices are intended to be integrated in a custom-made external lower limb prosthesis to ensure the function of the foot in patients with unilateral or bilateral lower limb amputation or deficiencies (transtibial/transfemoral amputation, knee/hip disarticulation, congenital limb deficiencies). They are intended for patients who would benefit from the dynamic toe response. This device is indicated for patients with moderate to high activity level (K3 to K4) for walking and physical activities without excessive overload.

RUSH ROVER maximum weight (load carrying included): 166 kg / 365 lbs

RUSH Kid maximum weight (load carrying included): 61 kg / 135 lbs

(See table §3)

6. CLINICAL BENEFITS

- Walking comfort
- Possible ambulation on uneven ground
- Ultra Low Profile

7. ACCESSORIES AND COMPATIBILITY

An appropriate foot shell must be mounted on the foot (refer to our catalog).

The foot includes a male pyramid connection designed to be compatible with standard female pyramid connectors (refer to our catalog).

8. ALIGNMENTS

Bench alignment


Recommendation: Add a 10 mm ($\frac{3}{8}$ ") wedge under back of the foot during bench alignment.

Sagittal plane:

- Define the appropriate socket flexion and heel height.
- Position the load line so that it falls through the center of the foot pyramid connection.

Frontal plane:

- Define the appropriate socket adduction/abduction.
- Position the load line so that it falls through the midline of the foot in a neutral mediolateral position.

Static alignment

Due to the shape of the rocker sole (bottom blade), patients may experience a new sensation when seeking the midpoint of the foot. The rocker sole allows patients to find their own comfortable static or standing position.

Adjustment of the anteroposterior set screws at the proximal connector is the more appropriate place to make changes regarding plantarflexion or dorsiflexion.

⚠ The heel wedges must not be introduced at this stage of alignment.

Dynamic alignment

The use of a slide adapter is recommended for dynamic alignment purposes, as it provides the best solution to troubleshooting the most common alignment issues listed below:

- Hard or soft heel
- Hard or soft toe
- Varus or valgus movements during stance phase

Once the optimal relative socket/foot alignment is determined, use plantarflexion or dorsiflexion at the proximal foot adapter to optimize heel-to-toe comfort and energy return.

9. ASSEMBLY

After dynamic alignment, tighten the pyramid adjustment screws according to the specifications of the connector manufacturer. Secure pyramid adjustment screws with a thread locking adhesive (e.g., Loctite 242).

Spectra sock

A Spectra sock is included to protect the foot shell and fiberglass components and minimize noise. It must be placed over the keel and sole plates before mounting the foot shell.

Foot shell

To install and remove the foot shell, use a shoehorn to prevent damage of the foot module.

- ⚠ Never remove the foot from the foot shell by pulling manually. Never use a screwdriver or any other inappropriate instrument to remove it. This could damage the foot.

10. TROUBLESHOOTING

- ⚠ If the patient notices any abnormal behavior or feel any changes in the characteristics of the device (noise, play, excessive wear...), or if the device has received a severe impact, they should stop using the device and consult their prosthetist.

11. WARNINGS

- ⚠ In case of damaged packaging, check the integrity of the device.
- ⚠ Never use the foot without a foot shell and a Spectra sock.
- ⚠ Never loosen the pyramid fastening screws.
- ⚠ The patient must inform their prosthetist if they gain or lose significant weight.
- ⚠ Make sure that the foot and inside of the foot shell are free of impurities (e.g., sand). The presence of impurities causes the fiberglass parts and the foot shell to wear out. Clean the foot according to the instructions (see §15).
- ⚠ After swimming or using in water, the foot including its foot shell must be cleaned (see §15).
- ⚠ Failure to follow the instructions for use is dangerous and will void the warranty.

12. CONTRAINDICATIONS

- ⚠ Do not use for a patient whose maximum weight (load carrying included) may exceed 166 kg / 365 lbs (RUSH ROVER) or 61 kg / 135 lbs (RUSH Kid).
- ⚠ Do not use for a patient who do not meet the requirements of an activity level of K3 or higher.
- ⚠ Do not use for activities associated with a risk of significant impact or excessive overloading.

13. SIDE EFFECTS

There are no known side effects directly associated with the devices.

Any serious incident that has occurred in relation to the devices should be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

14. MAINTENANCE AND CONTROL

No maintenance operation such as lubrication, work on the screws or other parts is required.

The foot module must be inspected by the prosthetist at least every six months. Inspections at shorter intervals are required if the user is more active.

The Spectra sock and the foot shell must be replaced by the prosthetist at regular intervals, depending on the patient's level of activity. If these parts are damaged, it can lead to premature foot wear.

15. CLEANING

Remove the foot shell and the Spectra sock, rinse the foot in clear water, clean it with neutral soap and dry carefully.

The foot shell can be cleaned with a damp cloth or sponge. It must be dry before reuse.

 The devices are not resistant to solvents. Exposure to solvents may cause damage.

16. ENVIRONMENTAL CONDITIONS

Temperature range for use and storage: -20°C to 60°C [-4 to 140° F]

Relative air humidity: no restrictions


Waterproof: the devices are resistant to fresh, sea, and chlorinated water.

 The foot shell is not resistant to ultraviolet light (UV). Do not store it under direct sunlight.

17. DISPOSAL

The devices are made of fiberglass composite material and metal. The devices and their packaging must be disposed of in accordance with local or national environmental regulations.

18. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of 1 st declaration
	Authorized representative in the European Union		Single patient, multiple use		

19. REGULATORY INFORMATION

These products are CE-marked medical devices and are certified as conforming with Regulation (EU) 2017/745.

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PROTEOR SAS